

INDIAN SCHOOL MUSCAT SENIOR SECTION DEPARTMENT OF CHEMISTRY CLASS XII CHEMISTRY PRACTICAL CONTENT BASED EXPERIMENTS

Experiment Number: 11

Date: -----

ANALYSIS OF FOOD

Aim: To identify the given samples as Protein and Carbohydrate.

	OBSERVATION FOR SAMPLE A & B		
EXPERIMENT			INFERENCE
	Α	B	
Treat the given sample solutions with few drops of Biuret solution .	Purple colouration	No characteristic observation	Sample A is protein
Acidify dilute solution of samples with con. HNO ₃ acid	Yellow precipitate.	No characteristic observation	Sample A is protein
Treat dilute solution of samples with few drops of Millon's reagent	White precipitate turns red on heating.	No characteristic observation	Sample A is protein
Mix 1ml of the given samples with 2ml of Fehling's solution (1 ml each of A & B) and heat on a water bath.	No characteristic observation	A reddish brown precipitate is formed.	Sample B is carbohydrate
Mix 1ml of the given samples with 1ml of Tollen's reagent * and heat on a water bath.	No characteristic observation	A silver mirror is formed on the inner walls of the test tube	Sample B is carbohydrate

* <u>Preparation of Tollen's Reagent</u>: Wash a test tube with a little NaOH solution and take about 1 ml of Silver nitrate (AgNO₃) solution. A black precipitate formed is just dissolves in minimum quantity of NH₄OH solution (Add drop wise carefully).

<u>Result:</u> The given sample A is Protein and sample B is ------



INDIAN SCHOOL MUSCAT SENIOR SECTION DEPARTMENT OF CHEMISTRY CLASS XII CHEMISTRY PRACTICAL CONTENT BASED EXPERIMENTS



Experiment Number: 12

ANALYSIS OF ORGANIC FUNCTIONAL GROUPS

Date: -----

Aim: To identify the functional group present in the given sample of organic compound.

	Experiment	Observation	Inference		
	TESTS FOR ALDEHYDES (-CHO group) (Use acetaldehyde)				
1.	Mix 1ml of the given sample with 2ml of Fehling's solution $(A + B)$ and heat on a water bath	A reddish brown precipitate is formed.	Presence of –CHO group.		
2.	Mix 1ml of the given sample with 1ml of Tollen's reagent and heat on a water bath.	A silver mirror is formed on the inner walls of the test tube.	Presence of –CHO group is confirmed.		
	TESTS FOR CARBOXYLIC ACID (-COOH) GROUP (Use acetic acid)				
1.	Add a little of a saturated solution of NaHCO ₃ to the given sample.	Brisk effervescence of a colourless odourless gas, which turns clear limewater milky.	Presence of carboxylic acid group.		
2.	Mix 1ml of the given sample with 1 ml of ethanol and 1 drop of con. H ₂ SO ₄ acid. Heat the mixture on a boiling water bath for about 5 minutes. Remove the test tube from the water bath, pour the contents into a beaker containing about 25ml water (tap water) and note the smell.	A pleasant fruity smell of ester is evolved.	Presence of carboxylic acid is confirmed.		
	TESTS FOR AMINO (-NH2) GROUP (Use aniline)				

1.	Shake 2 drops of the given organic compound with 2 ml of dil. HCl.	The compound dissolves.	Presence of –NH ₂ group.			
2.	Take 1 ml each of the given organic sample, con.HCl, NaNO _{2(aq)} and Alkaline β -naphthol solution in four different test tubes. Cool them below 5 ^o C in an ice bath. Mix them in the following order (Stir with a glass rod after each addition). First HCl + NaNO ₂ followed by Aniline and finally β -naphthol and stir for 5 minutes without removing the reaction test tube from ice bath.	A red orange coloured dye (Precipitate) is obtained.	Presence of –NH ₂ group is confirmed.			
	TESTS FOR PHENOLIC (- OH) GROUP (Use Phenol)					
1.	Treat the samples with neutral Ferric chloride solution (Add 10 ml of water)	Violet colouration	Presence of Phenol			
2.	Phthalic Test: - Take 2-3 flakes of phthalic anhydride and 1ml of the sample in a test tube. Add 2-3 drops of con. H_2SO_4 and heat on boiling water bath for about 2-3 minutes and pour the contents into dil.NaOH solution taken a beaker.	Pink colouration develops	Presence of Phenol			

Result

The given organic compound contains the functional group ------ ()
